quattro[®].mini

Leading in Orthodontic Implants

psm-medical.com

Quality is never an accident.

It is always the result of high intention, sincere effort, intelligent direction and skillful execution; it represents the wise choice of many alternatives.

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Introduction | Einführung

Orthodontic Mini-Anchor System

The use of Mini-Anchor Screws or Temporary Anchorage Devices (TAD's) has become a standard of care within the orthodontic community over the past few years. First offered in 2000 the **quattro**® TAD design has continuously been changed and improved to become the **quattro**® TAD family of today. The **quattro**® TAD's even after more than two decades of production are still unmatched in terms of versatility, quality and maximal success rate. Since year 2000 far more than 200.000 **quattro**® TAD's have been placed in more than 70 countries worldwide. With the **quattro**®.mini screw a new even smaller bracket head design is launched in 2012 offering a new dimension in patient comfort based on the outstanding **quattro**® design (see below QR Video and page 11). The **quattro**®.mini System is ideally enhanced with the unique **BENEFit**®/**BENEplate** System (developed in cooperation with University of Düsseldorf, Germany) offering exchangeable abutments and connecting plates. This opens up new possibilities in palatal based appliances popularly used for molar distalization and many more indications where safe skeletal anchorage is required.

Features:

- Clinically proven applications
- Biocompatible, high-strength Titanium
- Self-drilling & self-tapping technique
- Reliable & secure anchorage
- Innovative screwhead designs
- Stronger conical body design
- Deeper threads & denser pitch design
- Available in multiple sizes

Benefits:

- Provides improved higher strength
- Immediate load-bearing of heavy forces
- Easy to combine with applications of orthodontic accessories
- Optimal mechanical retention
- Avoids gingival impingement
- Reduced treatment time
- Easy to insert and remove
- Variable treatment spectrum

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quattro®RH micro

quattro®mini Aligner

RH: Round head design | Runder Kopf Design
BH: Bracket head design | Bracket Kopf Design

* Developed in cooperation with: | Entwickelt in Zusammenarbeit mit:
Dr. Nicola Derton Studio Odontoiatrico Specialistico, Conegliano (TV), Italy

One Screw – Many Options

The **quattro®mini** TAD system offers **4 different screwhead designs** which are perfectly adapted to the Orthodontist's appliances used daily. All **quattro®mini** screws (1.3 mm to 2.3 mm) feature the identical and well-proven thread design for easy insertion and most secure and satisfactory treatment results.

Based on the successful quattro Screw the **quattro®mini** screw has been reduced in size in the head and neck area to cope with the most difficult anatomical conditions. The design of the self-drilling thread is identical to the **PSM quattro®** and **BENEFIT®** screws which have demonstrated outstanding primary stability and extremely low failure rates for many years. The very sharp threads and the optimized thread pitch facilitate the insertion tremendously and make the **quattro®mini** the first choice for Standard Mini anchor screws.

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
Optimized thread design leads to easy insertion without force (cork screw principal).

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
Illustrations are for illustrative purposes only and are not to scale, except where scale is indicated.

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
33-18003

P1


Screwholder, with stop, DENTAL,
for quattro[®]mini

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The instruments can be used for several **PSM Medical** products. In order to find the right instruments for your use, the sterile screwholders are matching in colors with the correct drill.

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33-18244

P1

The extension of 16 mm works for all rotary instruments of the **quattro®** and **BENEFIT®** System with a latch type DENTAL attachment (max. 600 rpm).

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
Drill Stop | Bohrerstopp


The drill stop prevents the clinician from unintentionally drilling too deep. They assure an accurate working length from 7 to 13 mm in 2 mm increments. By sliding the stop sleeve over the appropriate **PSM** drill, the working length is adjusted to the desired depth. Removal of the sleeve is easy by using the buffer.

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
33-18240

Screwdriver | Schraubendreher

The manually turned unit is connected to your existing contra-angled handpiece (blue 1:1) that enables its use like an angled screwdriver.

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10-63025

P1

Manually turned unit mod. to Pauls, with adjustable torque from 0 – 40 Ncm

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

DENTAL attachment
DENTAL Ansatz

Battery Screwdriver – iSD900 | Akku Schraubendreher – iSD900

Cordless Battery Screwdriver with Torque Calibration System for the controlled fixation of PSM Mini screws.

The iSD 900 will guarantee precise torque control due to the unique – easy to use – torque calibration system (10 – 40 Ncm). Time-saving and precise, riskfree insertion and tightening of screws and abutments. The complete set contains a reliable and user friendly handpiece for Dental attachments. Reverse mode is indicated with sound. Standard rechargeable AAA Ni-Mh batteries can be used. → We recommend using the original batteries (Ni-MH rechargeable batteries).

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

PSM Digital Media Library – Your Real-Time Source for Orthodontic Excellence

Welcome to the official **PSM Medical Media Library**.

This central hub provides clinicians, distributors, partners and interested parties with instant, 24/7 self-service access to our complete portfolio.

- Always Up-to-Date: Our online documents (ePublications) are maintained in real-time.
- Verified Content: Every publication is officially approved by PSM Management, ensuring verified and revision-controlled information.
- Smart Access: No special software is needed. Access Catalogues, Guidelines, Certificates, Videos, Scientific Publications etc. on any device via a dedicated URL.
- Responsibility Note: We provide a PDF download option for offline use. However, please note that only the online version is part of our automated revision service. To ensure you are using the latest clinical data, we recommend referencing the live link.

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Notes | Notizen

Specifications and design are subject to change without prior notice. All pictures, photos and product descriptions are the property of PSM Medical GmbH. Utilization and copies by third parties have to be authorized by PSM Medical GmbH. All rights reserved! All trademarks and brands mentioned in this brochure are subject to unrestricted regulations of the trade mark laws valid and right of ownership of the registered owner.

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PSM North America, Inc.