



RMO® PurePak™ Wire Products

RMO's PurePak™ combines individual wire packaging with the highest quality selection of preformed arch wires. PurePak™ is available in all RMO® Thermalloy®, Thermalloy®Plus, Orthonal®, and Bendaloy® product lines, including Natural, Ideal, and Expanded arch forms.



Each arch wire is individually packaged and sealed in RMO's PurePak™ for enhanced product safety



PurePak™ protects against airborne contaminants and cross-contamination

Tear away tab enclosure for easy access to individually sealed wires

RMO's PurePak™ containers feature color-coded labels for easy identification of your wire inventories

Note: The PurePak™ system is not intended for autoclave or any other methods of sterilization.

RMO® Pre-Torqued Arch Wires



Size Inch	Maxillary 34mm	Maxillary 38mm	Mandibular 28mm
RECTANGULAR			
.016 x .022	A07803	A07806	A07800
.017 x .025	A07804	A07807	A07801
.019 x .025	A07805	A07808	A07802



Size Inch	Maxillary	Mandibular	PurePak™ Maxillary	PurePak™ Mandibular
ROUND				
.016	A07121	A07130	A07121-1	A07130-1
.018	A07122	A07131	A07122-1	A07131-1
SQUARE				
.016 x .016	A07124	A07133	A07124-1	A07133-1
RECTANGULAR				
.016 x .022	A07125	A07134	A07125-1	A07134-1
.017 x .025	A07126	A07135	A07126-1	A07135-1
.019 x .025	A07128	A07137	A07128-1	A07137-1
.021 x .025	A07129	A07138	A07129-1	A07138-1



Size Inch	Natural				Ideal			
	Maxillary	Mandibular	PurePak™ Maxillary	PurePak™ Mandibular	Maxillary	Mandibular	PurePak™ Maxillary	PurePak™ Mandibular
ROUND								
.013	A07310	A07330	A07310-1	A07330-1	A07360	A07380	A07360-1	A07380-1
.014	A07300	A07320	A07300-1	A07320-1	A07350	A07370	A07350-1	A07370-1
.015	A07311	A07331	A07311-1	A07331-1	A07361	A07381	A07361-1	A07381-1
.016	A07301	A07321	A07301-1	A07321-1	A07351	A07371	A07351-1	A07371-1
.017	A07312	A07332	A07312-1	A07332-1	A07362	A07382	A07362-1	A07382-1
.018	A07302	A07322	A07302-1	A07322-1	A07352	A07372	A07352-1	A07372-1
.020	A07313	A07333	N/A	N/A	N/A	N/A	N/A	N/A
SQUARE								
.016 x .016	A07303	A07323	A07303-1	A07323-1	A07353	A07373	A07353-1	A07373-1
.020 x .020	A07315	A07335	N/A	N/A	N/A	N/A	N/A	N/A
RECTANGULAR								
.014 x .025	A07314	A07334	N/A	N/A	N/A	N/A	N/A	N/A
.016 x .022	A07304	A07324	A07304-1	A07324-1	A07354	A07374	A07354-1	A07374-1
.017 x .025	A07305	A07325	A07305-1	A07325-1	A07355	A07375	A07355-1	A07375-1
.018 x .025	A07306	A07326	A07306-1	A07326-1	A07356	A07376	A07356-1	A07376-1
.019 x .025	A07307	A07327	A07307-1	A07327-1	A07357	A07377	A07357-1	A07377-1
.021 x .025	A07308	A07328	A07308-1	A07328-1	A07358	A07378	A07358-1	A07378-1



Size Inch	Natural				Ideal				Expanded	
	Maxillary	Mandibular	PurePak™ Maxillary	PurePak™ Mandibular	Maxillary	Mandibular	PurePak™ Maxillary	PurePak™ Mandibular	Universal Max/Man	Universal PurePak™
ROUND										
.013	A07700	A07720	A07700-1	A07720-1	A07740	A07760	A07740-1	A07760-1	E07010	E07010-1
.014	A07701	A07721	A07701-1	A07721-1	A07741	A07761	A07741-1	A07761-1	N/A	N/A
.015	A07702	A07722	A07702-1	A07722-1	A07742	A07762	A07742-1	A07762-1	E07011	E07011-1
.016	A07703	A07723	A07703-1	A07723-1	A07743	A07763	A07743-1	A07763-1	N/A	E07012-1
.018	A07704	A07724	A07704-1	A07724-1	A07744	A07764	A07744-1	A07764-1	N/A	E07013-1
.020	A07710	A07730	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A
SQUARE										
.016 x .016	A07705	A07725	A07705-1	A07725-1	A07745	A07765	A07745-1	A07765-1	N/A	E07019-1
.020 x .020	A07712	A07732	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A
RECTANGULAR										
.014 x .025	A07711	A07731	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A
.016 x .022	A07706	A07726	A07706-1	A07726-1	A07746	A07766	A07746-1	A07766-1	N/A	E07020-1
.017 x .025	A07707	A07727	A07707-1	A07727-1	A07747	A07767	A07747-1	A07767-1	N/A	E07021-1
.018 x .025	A07708	A07728	A07708-1	A07728-1	A07748	A07768	A07748-1	A07768-1	N/A	E07022-1
.019 x .025	A07709	A07729	A07709-1	A07729-1	A07749	A07769	A07749-1	A07769-1	N/A	N/A



Size Inch	Natural				Ideal				Expanded	
	Maxillary	Mandibular	PurePak™ Maxillary	PurePak™ Mandibular	Maxillary	Mandibular	PurePak™ Maxillary	PurePak™ Mandibular	Universal Max/Man	Universal PurePak™
ROUND										
.012	A07011	A07031	A07011-1	A07031-1	A07400	A07420	A07400-1	A07420-1	N/A	N/A
.013	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	E07000	E07000-1
.014	A07000	A07020	A07000-1	A07020-1	A07401	A07421	A07401-1	A07421-1	N/A	N/A
.015	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	E07001	E07001-1
.016	A07001	A07021	A07001-1	A07021-1	A07402	A07422	A07402-1	A07422-1	N/A	E07002-1
.018	A07002	A07022	A07002-1	A07022-1	A07403	A07423	A07403-1	A07423-1	N/A	E07003-1
SQUARE										
.016 x .016	A07004	A07024	A07004-1	A07024-1	A07405	A07425	A07405-1	A07425-1	N/A	E07004-1
.018 x .018	A07007	A07027	A07007-1	A07027-1	A07406	A07426	A07406-1	A07426-1	N/A	N/A
RECTANGULAR										
.016 x .022	A07005	A07025	A07005-1	A07025-1	A07407	A07427	A07407-1	A07427-1	N/A	E07005-1
.017 x .025	A07006	A07026	A07006-1	A07026-1	A07408	A07428	A07408-1	A07428-1	N/A	E07006-1
.018 x .025	A07008	A07028	A07008-1	A07028-1	A07409	A07429	A07409-1	A07429-1	N/A	E07007-1
.019 x .025	A07009	A07029	A07009-1	A07029-1	A07410	A07430	A07410-1	A07430-1	N/A	N/A
.021 x .025	A07010	A07030	A07010-1	A07030-1	N/A	N/A	N/A	N/A	N/A	N/A



THERMALLOY®

- Lower tooth moving forces
- Ideal for significant crowding / high irregularity cases
- Thermal-activated nickel-titanium
- Excellent resiliency and exhibits true thermal performance
- Soft at room temperature for easy ligation
- Outstanding shape integrity
- Gentle forces for greater patient comfort
- Very responsive to chilling

THERMALLOY+® PLUS

- Lower to moderate tooth moving forces
- Thermal-activated nickel-titanium
- Optimal midrange load characteristics
- Slightly thermal at room temperature for easy ligation
- Responsive to chilling
- Excellent resiliency

Orthonol®

- Consistent moderate forces
- Superelastic nickel-titanium
- High flexibility and resiliency
- Ideal for mid-stage treatment
- Responsive to chilling
- Less patient discomfort than higher force wires

bendaloy®

- Moderate to higher forces
- Nickel-free beta-titanium
- Eliminates nickel sensitivity concerns
- Solid bending characteristics and excellent formability
- Ideal for mid to finishing stages of treatment



PREFORMED ARCH WIRES
by Rocky Mountain® Orthodontics



The World's Oldest Synergistic, Bioprogressive®, Breathing Enhancement Orthodontic Company®.
To order, please contact your RMO® Sales Representative or call 800.525.6375

