



S-line™

S-line™ SL ceramic bracket


For the highest demands on esthetics, comfort and treatment efficiency



S-line™ is a passively self-ligating all-ceramic bracket manufactured using nanoceramic technology. Thanks to the high-purity fine-grained ceramic particles that are processed and refined in this process, a transparency is achieved that is unparalleled on the market. The resulting exceptionally high esthetics, coupled with an innovative patented design and the best features that a state-of-the-art SL ceramic bracket can offer today, make S-line™ a true masterpiece of appearance and functionality.

The advantages at a glance

- 30% more material strength and transparency due to manufacturing using nanoceramic technology and thanks to a special finishing process
- outstanding esthetics that harmonize perfectly with the natural tooth shade
- strongly rounded design and 10 % reduced profile height for highest wearing comfort
- patented, durable S.N.P. closure mechanism for reliable opening and closing at all times
- patented base adapted to the tooth curvature with inclined undercuts for optimal bonding and easy debonding
- minimized bow friction due to the rounded design of slot base and slot edges
- double wing design for maximum flexibility during treatment
- easy rotation control thanks to the wide and stable sliding clip



More information

MBT*, slot size .022"

	Maxillary	Art.-Nr.		torque	angulation	offset
		R	L			
●	Central	DL265111	DL265121	17°	4°	0°
●	Lateral	DL265112	DL265122	10°	8°	0°
●	Cuspid W/H	DL265113H	DL265123H	0°	8°	2°
●	1st Bicuspid W/H	DL265114H	DL265124H	-7°	0°	2°
●	2nd Bicuspid W/H	DL265115H	DL265125H	-7°	0°	2°

	Mandibular	Art.-Nr.		torque	angulation	offset
		R	L			
●	Anterior	DL265131	DL265141	-6°	0°	0°
●	Anterior	DL265131	DL265141	-6°	0°	0°
●	Cuspid	DL265133	DL265143	0°	3°	2°
●	1st Bicuspid	DL265134	DL265144	-12°	0°	0°
●	2nd Bicuspid	DL265135	DL265145	-17°	0°	0°

Set Max./Man. 5-5 (Upper 3-5 W/H) DL265171 *MBT is a registered trademark of 3M Unitek.

Roth, slot size .022"

	Maxillary	Art.-Nr.		torque	angulation	offset
		R	L			
●	Central	DL265411	DL265421	12°	5°	0°
●	Lateral	DL265412	DL265422	8°	6°	0°
●	Cuspid W/H	DL265413H	DL265423H	0°	11°	2°
●	1st Bicuspid W/H	DL265414H	DL265424H	-7°	0°	2°
●	2nd Bicuspid W/H	DL265415H	DL265425H	-7°	0°	2°

	Mandibular	Art.-Nr.		torque	angulation	offset
		R	L			
●	Anterior	DL265421	DL265431	-1°	0°	0°
●	Anterior	DL265422	DL265431	-1°	0°	0°
●	Cuspid	DL265423H	DL265433	-11°	3°	2°
●	1st Bicuspid	DL265424H	DL265434	-12°	0°	0°
●	2nd Bicuspid	DL265425H	DL265435	-17°	0°	0°

Set Max./Man. 5-5 (Upper 3-5 W/H) DL265471

